

# Report from diagnostics group

# **Nido Landen presentation: 6 x-ray diagnostic techniques for WDM**

**High energy flash x-ray radiography**

**K-alpha**

**Bremsstrahlung**

**EXAFS (extended x-ray absorption fine structure) measure photon energy spectrum to measure temperature (0.1 to 1 eV)**

**Diffraction (0.1 – 1 eV) measures dynamic response of crystal lattice**

**Coherent scattering - 2.4 eV ion structure**

**X-ray Thomson scattering (1-100 eV)**


**Collisionally broadened plasma resonance (1-300 eV)**

## Discussion on diagnostics

- **4- ways to measure temperature (R. More)**
  - Hydrodynamic release: x-rays, lasers, optical imaging
  - Electrical conductivity: need to worry about charged particles, contact potential. → Sandwich, perhaps Terahertz waves
  - Optical emission: optical laser probe/polarimetry; pyrometer has problem with thin optical skin depth  $\sim 100$  Angstrom
  - X-ray diffraction (previous slide)
- **dE/dx: Na/Li source; two-beam – measure energy of Li or H+ beam beyond target, ion beam radiography**
- **X-pinch used as an x-ray source (D. Hammer; U Nevada Reno)**
- **Measure pressure using shock waves: pulsed pump & cw probe lasers**
- **Beam spot: GSI – K-alpha; NDCX, HCX – gas cloud or scintillator**

# Conceptual design of a generic WDM target for reproducible manufacture (original sketch R. More).


- Target to be mounted on a remote positioner/ carousel target changer.
- Simple optical/electrical diagnostics as indicated
- Attention to
  - Dimension and position
  - Inductance
  - Optical transmission



## Target chamber

- **Larry Grisham: What is effect of 5-15 T solenoid magnetic field on target and diagnostics?**
- **Need to model solenoid with access holes for multiple target chamber ports: on axis, 90 degrees, and 45 degrees upstream and downstream.**
- **Target manufacturing cost?**
- **H. Yoneda:**
  - **Diagnostic mirror in center of annular beam**
  - **Suggests possibility of including part of this work in US-Japan collaboration.**

**Target chamber – require multiple ports in presence of strong solenoid final focus.**



M. Leitner

**Schematic drawing of a solenoid final focus system for HEDP studies on NDCX-II.**

## Initial diagnostic suite roughly prioritized

- Streak camera (purchase)
- Fast optical pyrometer (GSI?)
- VISAR (purchase)
- X-ray source (x-pinch?) as backlighter
- X-ray pinhole camera
- Optical spectrometer
- UV spectrometer
- X-ray spectrometer
- Short pulse pump laser